



# WeldSight™ 1.13.0 Software

July 16, 2024

## Software version compatibility

WeldSight Software	1.13.0
System Update	12.338.2
FOCUS PX	
OpenView SDK	1.1.98.0
QuickScan iX PA+	
NDT device API	6.7.0.4
OmniScan X3 series	
MXU & OmniPC	5.18

## Improvements and bug fixes

- Corrected alignment for Dual Linear Array™ (DLA) probes when using the merge option
- Merging files now enabled for the ODF file extension
- The shortcut "SHIFT + right-click" for the Measurement cursor can be set properly
- The refresh is now enabled for the Measurement cursor display
- Other minor improvements and bug fixes